


MyKnee Surgical Planning Report

LONG AXIS



CAUTION

This case is based on CT data:

REMOVE FROM THE BONE THE CARTILAGE AND SOFT TISSUES COVERING THE CUTTING BLOCK CONTACT AREAS.

All measurements shown are from the bone and do not include the thickness of the cartilage.

CASE CODE J_MOL_LTK_AL_05091967
SURGEON Alonso Lisón
SURGERY DATE 25 Feb 2014
SURGICAL APPROACH Medial
PRODUCT GMK Primary
CUTTING BLOCKS STD




PREOPERATIVE DATA [°]

HKA 177.0
 Femur Valgus (from bone) 3.0
 Tibia Varus (from bone) 3.5
 Tibia Posterior Slope 6.5
 Epic Axis vs Post Cond 2.5

LEFT TOTAL KNEE | DEFAULT | CHANGED



Femoral Implant Size 6 Narrow 5
 Tibial Implant Size 6

FEMUR



	DEFAULT	CHANGED	COMMENTS
FEMORAL RESECTIONS [mm]			
Lateral Posterior Cut	5.5	6.5	
Medial Posterior Cut	9.0	10.0	
Lateral Distal Cut	5.0		
Medial Distal Cut	8.0		
FEMORAL ANGLES [°]			
Valgus	0.0		
Extension	0.0		
ROTATION [°]			
External Rotation vs. post. cond.	4.0		

TIBIA



	DEFAULT	CHANGED
TIBIAL RESECTIONS [mm]		
Lateral Tibial Cut	10.0	
Medial Tibial Cut	6.5	
TIBIAL ANGLES [°]		
Valgus	0.0	
Posterior Slope	3.0	